













Φ3.3mm

Φ4.1mm

Φ4.8mm



SPΦ3.3NN

SΦ3.3RN

SPΦ3.3RN

TEΦ3.3RN

SΦ4.1RN

SPΦ4.1RN

TEΦ4.1RN

SΦ4.8RN

SPΦ4.8RN












































SΦ4.8WN

SPΦ4.8WN

TEΦ4.8WN

S-L-A Implants

6 mm														
8 mm	011010	013010	012010	014010	013070	012060	014050	013130	012110	013180	012160		014090	
10mm	011020	013020	012020	014020	013080	012070	014060	013140	012120	013190	012170		014100	
12mm	011030	013030	012030	014030	013090	012080	014070	013150	012130	013200	012180		014110	
14mm	011040	013040	012040	014040	013100	012090	014080	013160	012140					
16mm		013050			013110									



NN=Narrow Neck RN=Regular Neck WN=Wide Neck SP=Standard Plus SP=Standard TE=Tapered Effect S-L-A=Sand-blasted, Large grit, Acid-etched

NN(Narrow Neck)

Φ3.5mm



RN(Regular Neck)

Φ4.8mm



WN(Wide Neck)

Φ6.5mm



Octagonal system

Octagonal system






Solid system

Hybrid

Octagonal system

Solid system

Transfer parts



Prosthetic restoration

Cement-retained

Screw-retained

Screw-retained or cement-retained

Cement-retained

Cement-retained

Bar

Retentive anchor




Screw-retained

Cement-retained



Cement-retained

Cement-retained




Abutments titanium



TiBase Abutment/ Temporary Abutment



Copings plastic



Auxiliary parts screws

